

## PRODUCT INFORMATION

日本語データシート

**Product Name :** DynaCompetent Cells DH5 $\alpha$   
(Previous name: Competent Cell DH5 $\alpha$ )

**Code No. :** DS220

**Size :** 100  $\mu$ l  $\times$  10

**Competency :**  $> 5 \times 10^8$  cfu/ $\mu$ g (pUC19)

**Supplied product :** SOC medium, 1 ml  $\times$  10

*This product is for research use only*



### Description :

DynaCompetent Cells DH5 $\alpha$  is a high-efficiency chemically competent cell from *E. coli* DH5 $\alpha$  strain (one of the standard strains for molecular biology research) and suitable for a wide variety of cloning applications. The DH5 $\alpha$  cell has mutation of  $\phi 80lacZ\Delta M15$  and lacks *lacI<sup>q</sup>* gene, which allows blue-white color screening of transformants with X-gal (IPTG is not required).

### Genotype of *E. coli* strain DH5 $\alpha$ :

*supE44,  $\Delta lacU169(\phi 80lacZ\Delta M15)$ , hsdR17, recA1, endA1, gyrA96, thi-1, relA1*

### Quality Control :

Transformation was carried out according to the method described in this Product Information using supercoiled pUC19 plasmid. Transformants were plated on LB plates containing 50  $\mu$ g/ml ampicillin. The efficiency was confirmed to be greater than  $5 \times 10^8$  cfu/ $\mu$ g.

### Storage condition :

Stable at -80 $^{\circ}$ C with little or no loss in transformation efficiency for 12 months from the date of receipt. Competent Cells are sensitive to variation in temperature. Must be stored at -80 $^{\circ}$ C. Upon receipt, store the DynaCompetent Cells DH5 $\alpha$  in a freezer at -80 $^{\circ}$ C directly from a dry ice shipping box and store SOC medium at room temperature or at -80 $^{\circ}$ C.

### Handling of competent cells :

- Competent cells are sensitive to mechanical shock. Excessive mixing should be avoided.
- After thawing competent cells on ice, cells tend to lose transformation efficiency gradually. Transformation should be started immediately following thawing cells on ice.
- Use of refrozen competent cells is not recommended.

### Composition of SOC medium supplied :

20 g/L	tryptone
5 g/L	yeast extract
0.5 g/L	NaCl
0.186 g/L	KCl
2.4 g/L	MgSO <sub>4</sub> , anhydrous
4 g/L	glucose

## PRODUCT INFORMATION

---

### Transformation Procedure :

- Materials to be supplied by user
  - LB plates with antibiotic
  - Ice bucket with ice
  - 15 ml sterilized-polypropylene culture tubes
  - 42°C water bath
  - 37°C shaker
  - Sterile spreaders
  - 37°C incubator

If blue-white screening is required to select transformants,

- 20 mg/ml X-Gal in dimethylformamide (DMF)

### ● Transformation

1. Thaw one tube of competent cells on ice. One tube contains 100 µl of cells for each transformation.
2. Add DNA sample\* directly into the competent cells and mix by flicking the tube.

\* The volume of DNA sample should not exceed 5 % of that of competent cells (i.e. for 100 µl of competent cells, use  $\leq 5$  µl).

3. Incubate the tube on ice for 20 minutes.
4. Heat Shock the cells by placing the tube in 42°C water bath for 45 seconds. Do not mix or shake.
5. Remove the tube from the 42°C bath and place it on ice for 2 min.
6. Transfer the cells to a 15 ml sterilized culture tube containing 0.9 ml of SOC medium (pre-warmed at room temperature to 37°C). Culture the cells at 37°C for 1 hr in a shaker.
7. Spread an aliquot of the cells onto an LB agar plate containing appropriate antibiotic.

If blue-white color screening is required, spread 25 µl of 20 mg/ml X-Gal onto an LB agar plate and allow the reagent to absorb 30 minutes prior to inoculating cells. As DH5α does not have *lacI<sup>q</sup>*, IPTG is not required for blue-white screening.

8. Incubate the plate at 37°C overnight.

### Reference:

Sambrook, J. and Russell, D.W. (2001) Molecular Cloning: A Laboratory Manual, 3rd ed., Cold Spring Harbor Laboratory Press, Cold Spring Harbor, NY.

### Related Products:

DS225	DynaCompetent Cells Jet DH5α
DS210	DynaCompetent Cells JM109
DS218	DynaCompetent Cells JM109 for Electroporation
DS228	DynaCompetent Cells DH5α for Electroporation