

PRODUCT INFORMATION

日本語データシート

Product Name : DynaCompetent Cells Jet DH5 α , Large
(Previous name: Jet Competent Cell (DH5 α), Large)
Code No. : DS225L
Size : DS225 \times 5 (100 μ l of cell \times 50, Recovery Medium 1 ml \times 50)
Competency : $> 2 \times 10^8$ cfu/ μ g (pUC19)



This product is for research use only

Description :

In transformation procedure, DynaCompetent Cells Jet DH5 α requires neither heat shock nor culture after heat shock. Transformation of the DynaCompetent Cells Jet DH5 α can be completed within approximately 10 minutes and its efficiency is higher than 2×10^8 cfu/ μ g with the provided Recovery Medium*. The time-saving procedure is a great benefit for researchers and experimenters. The strain of the Jet Competent Cell, DH5 α , is one of the standard strains as competent cells in molecular biology applications. The DH5 α cell has mutation of $\phi 80lacZ\Delta M15$ and lacks *laqI^q* gene, which allows blue-white color screening of transformants with X-gal (IPTG is not required).

*: Recovery Medium is prepared based on SOC medium.

Genotype of *E. coli* strain DH5 α :

supE44, $\Delta lacU169(\phi 80lacZ\Delta M15), hsdR17, recA1, endA1, gyrA96, thi-1, relA1$

Quality Control :

The DynaCompetent Cells Jet DH5 α was tested for transformation efficiency using supercoiled pUC19 plasmid according to the Jet Transformation Protocol described in this Product Information (LB plates containing 50 μ g/ml ampicillin) and its efficiency was confirmed to be greater than 2×10^8 cfu/ μ g.

Storage condition :

Stable at -80°C with little or no loss in transformation efficiency for 12 months from the date of receipt. Competent cells are sensitive to variation in temperature. Must be stored at -80°C. Upon receipt, store the DynaCompetent Cells Jet DH5 α in a freezer at -80°C directly from a dry ice shipping box and store Recovery Medium at 4°C (or -80°C). To avoid precipitation in Recovery Medium, slow freezing or freeze-thaw cycle (for example, storage around -20°C) should not be done, although transformation efficiency is not very affected by the precipitate.

Handling of competent cells :

- Competent cells are sensitive to mechanical shock. Excessive mixing should be avoided.
- After thawing competent cells on ice, cells tend to lose transformation efficiency gradually. Transformation should be started immediately following thawing cells on ice.
- Use of refrozen competent cells is not recommended.

PRODUCT INFORMATION

Transformation Procedure :

- Materials to be supplied by user

- LB plates with antibiotic
- Ice bucket with ice
- Sterile 1.5 ml tubes
- Sterile spreaders
- 37°C incubator

If blue-white screening is required to select transformants,

- 20 mg/ml X-Gal in dimethylformamide (DMF)

- Jet Transformation Protocol

1. Thaw one tube of competent cells on ice. One tube contains 100 µl of cells for each transformation.

2. Add DNA sample* directly into the competent cells and mix by flicking the tube.

* The volume of DNA sample should not exceed 5 % of that of competent cells (i.e. for 100 µl of competent cells, use ≤ 5 µl).

3. Incubate the tube on ice for 5 minutes.

4. Transfer the cells to a new 1.5 ml sterile tube containing 0.9 ml of Recovery Medium (pre-warmed at room temperature to 37°C), mix the tube contents by vortex for one second, and incubate the tube at room temperature for 5 minutes.

5. Spread all or an aliquot of the cells to an LB agar plate containing appropriate antibiotic.

If blue-white color screening is required, spread 25 µl of 20 mg/ml X-Gal on an LB agar plates and allow the reagent to absorb 30 minutes prior to inoculating cells. As DH5α does not have *lacI^q*, IPTG is not required for blue-white screening.

Note: It is especially important to absorb these solutions prior to inoculating cells for kanamycin or tetracycline selection. Do not mix cells with solutions of these reagents before inoculating to a plate.

6. Incubate the plate at 37°C overnight.

Reference:

Sambrook, J. and Russell, D.W. (2001) Molecular Cloning: A Laboratory Manual, 3rd ed., Cold Spring Harbor Laboratory Press, Cold Spring Harbor, NY.

Related Products:

DS230	DynaCompetent Cells JetGiga DH5α
DS220	DynaCompetent Cells DH5α
DS228	DynaCompetent Cells DH5α for Electroporation