

## PRODUCT INFORMATION

日本語データシート



- Product Name:** AllView PAGE Buffer®
- Code No:** DS520
- Size:** 500 ml (20 × stock solution)
- Storage:** Store at Room Temperature (Protect from light)
- Stability:** 12 months at Room Temperature

### Description

AllView PAGE Buffer® is a new type of running buffer for SDS-PAGE electrophoresis. This buffer has a remarkable feature that enable us to separate proteins with wide range of molecular weights in the basic Laemmli gel (Tris-HCl) similar to using “gradient gel”. The recommended acrylamide concentration is 6% for resolving gel and 3% for stacking gel to separate the entire range of about 10 kDa to 250 kDa. After electrophoresis, the polyacrylamide gel can be used directly for CBB staining, silver staining or western blotting.

### Protocol

1. Dilute the AllView PAGE Buffer® 20 times with ultrapure water.  
(e.g. Add 25 ml of AllView PAGE Buffer® to 475 ml of ultrapure water.)
2. Set the gel in electrophoresis chambers.
3. Fill the chambers with the 1 × AllView PAGE Buffer®.
4. Load your samples and molecular weight markers.
5. Start electrophoresis.

Gel	Voltage	Time
Mini gel (8 × 10 cm, 1 mm thick)	250 V (Constant voltage)	13~15 min

### Note:

- AllView PAGE Buffer® is suitable for Laemmli gels (see below “Recommended usage”).
- If precast gels are used, the optimal acrylamide concentration may be different.
- AllView PAGE Buffer® is not reusable.
- AllView PAGE Buffer® is 20 × stock solution. Please dilute to 1 × solution before use.
- If precipitation of SDS is observed, completely dissolve it in a water bath (about 37°C) before use.
- Optimal electrophoretic time is depending on acrylamide concentration, gel size and voltage.

Accordingly, it is recommended to monitor electrophoresis using a prestained MW Marker (e.g. DynaMarker® Protein MultiColor Stable II, Code#DM660).

- Sometimes the gel become too hot during electrophoresis with two gels at the same time or with a large size gel. To avoid heat generation, use the chilled AllView PAGE Buffer® and electrophorese in a cold room.

## PRODUCT INFORMATION


### Recommended usage

The recommended acrylamide concentration is 6% for resolving gel and 3% for stacking gel to separate the entire range of about 10 kDa to 250 kDa. Especially if the separation of low molecular weight proteins (10 kDa~30 kDa) is required, 10% resolving gel is recommended.

#### 1. Gel preparation (Laemmli's method)

**Table 1.** Recipes for polyacrylamide resolving and stacking gel (2 mini gels)

	Stacking gel (6 ml)	Resolving gel (15 ml)	
Gel percentage	3%	6%	10%
Ultrapure water	3.9 ml	8.05	6.05
1.5M Tris-HCl (pH8.8)		3.75	3.75
0.5M Tris-HCl (pH6.8)	1.5		
30% Acrylamide/Bis solution	0.6	3.0	5.0
10% SDS	0.06	0.15	0.15
TEMED	0.003 (3 $\mu$ l)	0.00375 (3.75 $\mu$ l)	0.00375(3.75 $\mu$ l)
10% APS	0.06	0.05	0.05




\* DynaMarker® Protein MultiColor Stable II, Code#DM660

## PRODUCT INFORMATION

### 2. Electrophoretic image

The mobility of proteins using AllView PAGE Buffer<sup>®</sup> is different from Laemmli electrophoresis running buffer (Tris-Glycine-SDS). By using AllView PAGE Buffer<sup>®</sup> and the Laemmli gels (Tris-HCl), a wide range of protein sizes can be separated like a gradient gels.



\* DynaMarker<sup>®</sup> Protein MultiColor Stable II, Code#DM660

### Related Products

Code	Product Name	Description
DM660	DynaMarker <sup>®</sup> Protein MultiColor Stable II	The DynaMarker <sup>®</sup> Protein MultiColor Stable II is a pre-stained protein molecular weight marker. The marker has a remarkable feature that it is possible to store at 4 °C.
DS500	QuickBlue Staining Solution	QuickBlue Staining Solution stains proteins in polyacrylamide gel after SDS gel electrophoresis (detection limit $\geq$ 8 ng of protein) with *Coomassie <sup>®</sup> -G250. All processes including washing and destaining processes can be performed in 1.5 hr.